


FRIABILITY TESTER GA - FT 2D



The GA-FT-2D is designed for Faster and Accurate operation. Operator can test according to USP/EP test protocols. Single or Dual drums can be used.

Supplied with complete validation documents (IQ/OQ). Weight of the tablets can be transferred directly from the balance before and after test easy calculation of the results.



GALAXY
INTERNATIONAL
+91 40 3580 7244

WWW.GICORP.IN INFO@GICORP.IN
PLOT NO 162, SRI RAMANPURAM STREET NO. 3
UPPAL, HYDERABAD, TELANGANA, INDIA, 500039

SYSTEM SPECIFICATIONS

- **No of Drums:**
 - 1 or 2 Drums
- **Test Methods:**
 - USP 1216 / EP 2.9.7 or any other method
- **Rotation Speed:**
 - 20 rpm - 50 rpm
- **Chamber:**
 - Made of SS304
- **Balance Connectivity:**
 - Available
- **Qualification:**
 - IQ-OQ and PV documentation available
- **MOC:**
 - Available in SS 304 & MS-PC

ADDITIONAL FEATURES

- External 3” dot Matrix printer connectivity (Epson/TVS)

